

Transducers

C5-2Q



2 - 5 MHz
Curved Array

P5-1Q



1 - 5 MHz
Phased Array

L12-5Q



5 - 12 MHz
Linear Array

E8-4Q



4 - 8 MHz
Endocavity
Curved Array

C5-2MQ



2 - 5 MHz
3D/4D Mechanical
Curved Array

P7-3Q



3 - 7 MHz
Phased Array

L17-7HQ



7 - 17 MHz
High Frequency
Linear Array

E10-3HQ



3 - 10 MHz
Endocavity Curved Array (200°)

MC8-4Q



4 - 8 MHz
Micro Convex Array

MC9-3TQ



3 - 9 MHz
Micro Convex Array

L17-7SQ



7 - 17 MHz
High Frequency Compact
Linear Array

E10-3BQ



3 - 10 MHz
Endocavity Curved Array (200°)

General Imaging Ultrasound in a Whole New Light.

Acclarix™ AX8
Diagnostic Ultrasound System



GESTURE-CONTROL USER INTERFACE

APPLICATION-BASED WORKFLOW

USER-FOCUSED DESIGN

UNPARALLELED VALUE

DEFINITIVE IMAGE QUALITY

Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com


© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.




Stunning Clarity

Sleek design and user-friendly features are just part of the story. The Acclarix AX8 is also designed to deliver amazing clarity and stunning image quality quickly—eliminating the need to tweak imaging controls.

- High fidelity, high-channel count architecture results in superb detail resolution, particularly at depth
- Tissue Adaptive Imaging (TAI) continuously and automatically optimizes imaging allowing more focus on the patient
 - In B-mode, TAI fine tunes multiple parameters to provide the best possible image quality
 - In Doppler, TAI automatically adjusts for flow state providing improved continuity, border detection and fill-in



L12-5Q
eVol.Flow



L12-5Q
Elastography




C5-2Q
Kidney



C5-2MQ
3D/4D



E10-3HQ
Uterus



P5-1Q
TDI

It's a Tour de Force

Redefining innovation through value and performance

It's the dawning of a new day in the world of compact ultrasound—and it turns out you really can have it all. The remarkable Acclarix AX8 Compact Ultrasound System delivers a powerhouse combination of features in a remarkably small package—making it the ideal choice for tight environments. The Acclarix AX8 has been designed from the ground up with a relentless focus on delivering unexpected levels of innovation and performance at a price point that is equally surprising.

- Definitive Imaging
- Distinctive Design
- Intelligent Workflow
- Intrinsic Quality

The Virtue of Value

Born of a vision to deliver meaningful design innovations that benefit the user, the Acclarix AX8 features a host of design breakthroughs that make day-to-day operation easy, fast and intuitive. The result is an elegant simplicity where form and function meet at the tips of your fingers. The Acclarix AX8 delivers unmatched value and performance across a broad range of applications.

- Abdomen
- Breast
- Prostate
- Difficult-to-Image
- Nerve
- Vascular
- OB/GYN
- Musculoskeletal
- Cardiac Screening
- Small Parts
- Intraoperative

Distinctive Design



Sleek compact design
for portable exam environments



Tilt-and-swivel monitor
for optimal viewing
during procedures



Dual touch screens
streamline exam workflow



Customizable control panel
to personalize
your user interface