


Streamlined Data Connectivity



Tailored for Advanced Cardiac Practice

SE-1202 Electrocardiograph



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Global Headquarters:
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen 518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:
EDAN Diagnostics, Inc. | 9918 Via Pazar, San Diego, CA 92126 | +1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

CE 0123
ENG-ECG-SE-1202
V1.0-20210519



Comprehensive Solutions for Cardiac Practice

SE-1202, as a next-generation ECG workstation, provides a comprehensive set of cardiac diagnosis tools, delivering cutting-edge technologies to support your clinical decisions in cardiac practice.



10.1" Touch UI



Internal Storage



Water-proof Keys



Pacemaker Detection

- 12-lead Standard Resting ECG
- 9-lead Resting ECG Tailored for Neonates
- Pharma Study
- Vector ECG & Timed Vector ECG (VCG&TVCG)

- Signal-average ECG (SAECG)
- Heart Rate Variability (HRV)
- ST-view (Intuitive Display of ST-segment)

The desirable 10.1" capacitive touch screen presents a gorgeous trace view, coming along with easy operations on an on-screen toolbar with customizable hotkeys.



Foldable screen with living hinge design.



The waterproof capacitive keys with multi-color breathing lights enable easy cleaning and sterilization, introducing a smooth input experience to your daily practice.



Exceptional Diagnosis Supports in Unexceptional Cases

Long-time Sampling

For certain cases in which a prolonged practice is required, an 1800 seconds recording could be applied, to help with either pretreatment diagnosis or post-treatment evaluations.

Auto Arrhythmia Detection

Arrhythmia trace segments could be highlighted on display with automatic saving strips for further analysis.

On-screen Diagnosis

Designed as an "All-in-One" workstation, the SE-1202 allows a complete paperless workflow with all measurements and diagnosis on screen, correctly and coherently. Any waveform segments can be amplified and measured. Any report can be edited and confirmed. Any two reports can be selected and compared.