

GW1 Gateway

Integration Solution




EDAN

A world of potential

GW1 Gateway Integration Solution

As a versatile informatics component, GW1 gateway software enables bi-directional communication between EDAN monitoring products and hospital information network via IHE standardized HL7 protocol or secondary development. It delivers seamless network access and data integration services for up to 256 EDAN patient monitors.



Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.