

Always Attend to
your Blood Pressure

SA-10 Ambulatory Blood
Pressure Monitor



SA-10 Ambulatory Blood Pressure Monitor

The Edan ABPM SA-10 is a 24-hour ambulatory blood pressure monitor with stylish design, intelligent software, state-of-the-art blood pressure monitoring technologies, developed to avoid white coat hypertension, masked hypertension and sustained hypertension. With the help of SA-10's powerful analysis software, clinicians will manage patients' cardiovascular risks with ease and make clinical decisions with confidence.



Self-locking structure



300NIBP
Memory



IP22*



Wireless transmission



<160g
Lightweight design



Glass panel

Feature Highlights

- Lightweight, compact and quiet for better patients compliance
- Large blood pressure measurement range from 10mmHg to 290mmHg
- Proprietary algorithm with excellent motion resistance performance
- Data protection from power outage
- Auto re-measurement on failed measurement attempt


*IP22 ingress protection rating is only with the carrying pouch together

Powerful Software for Comprehensive Analysis


With the aid of the SA-10's powerful analysis software, the clinicians can have a better insight into the variation of patient's blood pressure within 24 hours, easily identifying the periods when the blood pressure values are above limits. The informative statistical table provides a clear presentation of the measurement summary throughout the day.



Pie Chart Analysis



Circadian Analysis



Statistical Analysis



SA-10

Ambulatory Blood Pressure Monitor



A world of potential

Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-ECG-SA-10
V1.1-20240515